

LENTIS® | FEMTIS®
ACUNEX®



Viscoject Bio 2.2

Accuject 2.1-BL

LP604340C

LP604595

			Viscoject Bio 2.2	Accuject 2.1-BL
			LP604340C	LP604595
MONOFOCAL	LENTIS® L-302-1	0.0D to +35.0D	YES	NO
	LENTIS® L-312	-10.0D to +35.0D	YES	NO
	LENTIS® LS-312Y	0.0D to +35.0D	YES	NO
	LENTIS® L-303	0.0D to +35.0D	YES	NO
	LENTIS® L-313	-10.0D to +35.0D	YES	NO
	LENTIS® LS-313Y	0.0D to +35.0D	YES	NO
	LENTIS® L-323	+10.0D to +35.0D	YES	NO
	myLENTIS® LU-323	0.0D to +35.0D	YES	NO
	LENTIS® LU-814VR	0.0D to +35.0D	NO	YES
	FEMTIS® FB-313	+15.0D to +30.0D	YES	NO
	ACUNEX® AN6	+10.0D to +30.0D	NO	YES
Enhanced Mono.	LENTIS® QUANTUM L-333	+10.0D to +30.0D	YES	NO
	ACUNEX® QUANTUM AN6Q	+10.0D to +30.0D	NO	YES
TORIC	LENTIS® T _{plus} LS-313 T0-T6	+10.0D to +30.0D	YES	NO
	LENTIS® T _{plus} X LU-313 T TY	-10.0D to +35.0D	YES	NO
	myLENTIS® LU-323 T TY	0.0D to +35.0D	YES	NO
	LENTIS® T _{plus} LU-814 T	-10.0D to +35.0D	NO	YES
EDOF	LENTIS® Comfort LS-313 MF15	-10.0D to +36.0D	YES	NO
	FEMTIS® Comfort FB-313 MF15	+15.0D to +30.0D	YES	NO
	LENTIS® Comfort LU-814 MF15	0.0D to +36.0D	NO	YES
	ACUNEX® VARIO AN6V	+10.0D to +30.0D	NO	YES
EDOF TORIC	LENTIS® Comfort ^{toric} LS-313 MF15T0-T6	+10.0D to +30.0D	YES	NO
	FEMTIS® Comfort ^{toric} FB-313 MF15T0-T3	+15.0D to +30.0D	YES	NO
	ACUNEX® VARIO ^{TORIC} AN6VTO-T3	+10.0D to +30.0D	NO	YES
MULTIFOCAL	LENTIS® M _{plus} LS-313 MF20	-10.0D to +36.0D	YES	NO
	LENTIS® M _{plus} LS-313 MF30	-10.0D to +36.0D	YES	NO
	LENTIS® M _{plus} X LS-313 MF30	-10.0D to +36.0D	YES	NO
	LENTIS® M _{plus} X LS-313 MFX30	-10.0D to +36.0D	YES	NO
	LENTIS® M _{plus} LU-814 MF20	0.0D to +36.0D	NO	YES
	LENTIS® M _{plus} LU-814 MF30	0.0D to +36.0D	NO	YES
	FEMTIS® M _{plus} FB-313 MF30	+15.0D to +30.0D	YES	NO
	ACUNEX® VARIO ^{MAX} AN6VM	+10.0D to +30.0D	NO	YES
MULTIFOCAL TORIC	LENTIS® M _{plus} ^{toric} LU-313 MF15T TY	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} ^{toric} LU-313 MF20T TY	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} ^{toric} LS-313 MF30T0-T6	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} ^{toric} LU-313 MF30T TY	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} X ^{toric} LU-313 MFX30T TY	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} X ^{toric} LU-313 MF30T TY	0.0D to +36.0D	YES	NO
	LENTIS® M _{plus} ^{toric} LU-814 MF15T TY	0.0D to +36.0D	NO	YES
	LENTIS® M _{plus} ^{toric} LU-814 MF30T TY	0.0D to +36.0D	NO	YES
	FEMTIS® M _{plus} ^{toric} FB-313 MF30T0-T3	+15.0D to +30.0D	YES	NO
	ACUNEX® VARIO ^{MAX} ^{TORIC} AN6VMT0-T3	+10.0D to +25.0D	NO	YES

Recommended combination of IOL and Injector: YES NO